

Results of the Study – *A. brasiliensis* 9642

Test Microorganism	Carrier Type	Run Type	Replicate	Replicate CFU/carrier	Average CFU/carrier	Average Percent Reduction vs. Time Zero	Average Log ₁₀ Reduction vs. Time Zero
<i>A. brasiliensis</i> ATCC 9642	Glass	Control	1	4.00E+05	4.00E+05	N/A	
			2	4.00E+05			
			3	4.00E+05			
		Test	1	2.00E+01	7.00E+01	99.98%	3.76
			2	1.30E+02			
			3	6.00E+01			
	Vinyl	Control	1	7.00E+05	6.00E+05	N/A	
			2	7.00E+05			
			3	4.00E+05			
		Test	1	2.00E+01	9.67E+01	99.98%	3.79
			2	3.00E+01			
			3	2.40E+02			

