


Results of the Study – *S. aureus* 33592 (MRSA)

Test Microorganism	Carrier Type	Run Type	Replicate	Replicate CFU/carrier	Average CFU/carrier	Average Percent Reduction vs. Time Zero	Average Log ₁₀ Reduction vs. Time Zero
<i>S. aureus</i> ATCC 33592 (MRSA)	Glass	Control	1	6.10E+06	5.67E+06	N/A	
			2	6.20E+06			
			3	4.70E+06			
		Test	1	<1.00E+01	<1.00E+01	>99.9998%	>5.75
			2	<1.00E+01			
			3	<1.00E+01			
	Vinyl	Control	1	2.80E+06	1.14E+07	N/A	
			2	2.45E+07			
			3	6.90E+06			
		Test	1	<1.00E+01	<1.00E+01	>99.9997%	>5.51
			2	<1.00E+01			
			3	<1.00E+01			



*Note: the limit of detection for this assay is 1.00E+01 CFU/carrier. Values observed below the limit of detection are reported as <1.00E+01 in the table and as zero in the graph