


Results of the Study – *C. difficile* 43598 (Endospores)

Test Microorganism	Carrier Type	Run Type	Replicate	Replicate CFU/carrier	Average CFU/carrier	Average Percent Reduction vs. Time Zero	Average Log ₁₀ Reduction vs. Time Zero
<i>C. difficile</i> ATCC 43598	Glass	Control	1	1.35E+06	1.23E+06	N/A	
			2	1.13E+06			
			3	1.20E+06			
		Test	1	<1.00E+01	<1.00E+01	>99.9992%	>5.09
			2	<1.00E+01			
			3	<1.00E+01			
	Vinyl	Control	1	1.33E+06	2.90E+06	N/A	
			2	9.60E+05			
			3	6.40E+06			
		Test	1	1.50E+02	2.40E+02	99.992%	4.08
			2	6.00E+01			
			3	5.10E+02			



*Note: the limit of detection for this assay is 1.00E+01 CFU/carrier. Values observed below the limit of detection are reported as <1.00E+01 in the table and as zero in the graph